



# Duo System Premium and Duo System



## Repositionable Localization Systems


Localization wires designed for precise marking of breast lesions featuring a double hook wire design for excellent anchoring in tissue.


The Duo System Premium and Duo System are correctable localization systems for mammography. They offer doctors and patients top quality and reliability for the preoperative marking of suspicious tissue.

The localization wires come preloaded in 20-gauge thin, sharp puncture cannulas to ensure minimally invasive treatment. All localization wires are made of the biocompatible implant material, nitinol, and are available with different cannula lengths.

Features:

- Duo System Premium offers single-handed operation due to ergonomic handle that allows users to make markings and guide ultrasound transducer simultaneously.
- Flexible marking wire for intraoperative localization of findings.
- Secure anchoring in the tissue through double hooks with sharp tips.





### Ordering Information

Ref. no.	Description	Gauge	Diameter	Length	Qty/box
271770	Duo System Premium	19,5	0.95 mm	5 cm	10
271771	Duo System Premium	19,5	0.95 mm	9 cm	10
271773	Duo System Premium	19,5	0.95 mm	12 cm	10
271440	Duo System	20	0.95 mm	5 cm	10
271441	Duo System	20	0.95 mm	9 cm	10
271443	Duo System	20	0.95 mm	12 cm	10

MM 02-7000-03 ENG Rev. 03.

#### Contact Details, Nordics

customer@mermaidmedical.com  
mermaidmedical.com

#### Contact Details, NL

info@mermaidmedical.nl  
mermaidmedical.com

#### Contact Details, UK

customer@mermaidmedical.co.uk  
mermaidmedical.com