



HSG Catheter



HSG Catheters from CooperSurgical® are designed to diagnose blockages or abnormalities in the uterine cavity and fallopian tubes. They can be used for both hysterosalpingography (HSG) and hysterosonography.

These thin, flexible tubes are easily inserted into the cervix, where the latex-free concentric balloon helps secure the catheter in place. The catheters are compatible with both oil-based contrast dyes and saline-based solutions, making them suitable for both X-ray and ultrasound imaging.

Features:

- Flexible radiopaque polyurethane catheter shaft ensures visibility on X-ray for precise placement.
- Soft, latex-free concentric balloon placed 3 mm from the catheter tip to minimize patient discomfort.
- Atraumatic smooth round catheter tip allows easy guide wire introduction.
- Pre-set syringe with stopcock included for convenient balloon inflation.
- Available in two diameters: 5F and 7F.

Ordering Information

| Ref. no. | Description | Accepts Guidewire | Qty/box |
|----------|---------------------------------|-------------------|---------|
| 61-5005 | HSG Catheter 5F (Side delivery) | 0,018" | 10 |
| 61-5007 | HSG Catheter 7F (Side delivery) | 0,018" | 10 |

MM 02-7000-05 ENG Rev. 03

Contact Details, Nordics

customer@mermaidmedical.com
mermaidmedical.com

Contact Details, NL

info@mermaidmedical.nl
mermaidmedical.com

Contact Details, UK

customer@mermaidmedical.co.uk
mermaidmedical.com